

# Fluoroscopy



# • Definition

- Use of x-rays to produce a dynamic sequence of images ,which are viewed in real time




## Types:

- Direct viewing fluoroscopy
- Image intensifier
- Flat panel systems



# 1) Direct viewing fluoroscopy



- Oldest fluoroscopy
- Radiology directly viewed using zinc sulphide fluorescent screen
- screen backed with lead glass for protection
- Limitation: dim image
  - Must be carried in dark room, with very low contrast images (rods =scotopic vision)




Dr. Youssef Gamal

Fluorescent  
screen  
[www.ScoreTraining.net](http://www.ScoreTraining.net)

## 2) Image intensifier



# The Image Intensifier




- Components of image intensifier:

**1- Glass or ceramic envelope surrounded by metal housing:**

Function :

- A) Prevent light from getting in the tube
- B) Shield the device from magnetic field



## 2- Input screen :

- Consists of:




### A) Input phosphor layer:

- Laid down in a thin metal layer
- At the outer (X-ray beam) side of the input screen
- Always: cesium iodide
- Advantages:

- ☐ High absorption efficiency (About 60% of the X-ray photons in 0.1-0.4 thick phosphor layer)

Because k-edge of cesium = 36 kev , and I = 33 kev (just below the used effective energy)

- ☐ Low unsharpness with thicker phosphors (why?)




# Input screen (cont.)

- Size: 150 – 400 mm in diameter
- Smallest for fracture fixation
- Largest for angio




### 3) photocathode:




- Thin coating material on the inner side of the input screen
- It emits electrons when irradiated by light
- Number of electrons emitted  $\propto$  with intensity of light  $\propto$  with intensity of X-ray
- Commonly : antimony caesium
- It is maintained with high –ve voltage with respect to anode e.g 25kv

#### **4- electron focusing electrodes:**


- Metal rings within the tube
- Has + ve voltage with respect to photocathode
- Constrain electrons to travel directly to output screen so that: pattern of electron intensities falling on the screen are identical as (but minified) pattern of intensities produced in the input screen
- i.e. functions as electron lens
- If absent : electrons reach the output screen in disorganized pattern → no image formed





## 5- output screen:

- 25-35 mm in diameter (10 times smaller than input screen)




- Consists of:

**A) Thin phosphor layer (zinc cadmium sulphide  $\text{ZnCdS}$ ):**


- Converts pattern of electron intensities into light
- It must be so thin (few  $\mu\text{ms}$ ) to decrease light spread in the screen
- It Emits green light
- Brightness displayed in the output screen is  $\propto$  intensity of X-ray (X-ray doubled  $\rightarrow$  screen Brightness doubled) i.e.  $\gamma=1$

**B) Extremely thin layer of aluminum (0.5) covering the output screen**

- Function:
  - 1) Prevent light emitted backwards towards the tube from reaching photocathode causing cascade of electron emission that would completely whiteout the image
  - 2) Acts as anode with +ve voltage with respect to cathode e.g 25kv  $\rightarrow$  acceleration of electrons towards the output screen




# Summary of image intensifier



# Fluoroscopy image characteristics

## ***Brightness Gain:***

- Definition: Ratio of brightness of the output phosphor to that of the input phosphor = degree of intensification



- Brightness gain is Caused by two factors:

1) Flux gain :


Single light photon from the input phosphor → single electron from the photocathode → accelerated at 25Kev → many (about 50) light photons from output phosphor

2) Minification gain:

Intensification caused by decrease of the image size from input to output screen

Equal to ratio of the **area** of both screens =  $300^2/30^2 = 100$

**Overall brightness gain = Flux gain x Minification gain = 50 x 100 = about 5000**







# ***Conversion factor = Gx:***

$$\text{Conversion factor} = \frac{\text{brightness (luminance) at output phosphor (candela.m}^{-2}\text{)}}{\text{dose rate at the input surface of image intensifier (}\mu\text{Gy .s}^{-1}\text{)}}$$


- Used to describe the performance of image intensifier (as the gain is not measurable)
- Typical value = 25-30 candela.m<sup>-2</sup> (μGy .s<sup>-1</sup>)<sup>-1</sup>
- Gain and Gx decrease with equipment usage due to loss of detection efficiency of phosphors

- ***Magnification:***
- Image in the output phosphor is a minified version of the image
- To magnify the image:
  - Voltage of the intermediate electrodes is changed → electron focus moves backwards → central part of the image only fill the output phosphor
- Advantages: better resolution






Normal Operation




Magnification Mode

- Disadvantages: ↓ minification gain → ↓ image brightness
- To restore brightness → ↑ exposure factors with :  
↑ entrance skin dose
- All the image intensifiers offer at least one magnified fields of view




# Closed Circuit TV Systems

- Used to display the image of the output phosphor
- Consists of :
  - TV Camera
  - Camera control unit
  - Monitor
- Image intensifiers and television camera tubes are manufactured so that the output phosphor of the image-intensifier tube is the same diameter as the window of the television camera tube, usually 2.5 or 5 cm





# Image Intensifier – TV camera Coupling

- 1) Lens coupling
- 2) Fiber Optic coupling





# Lens Coupling

- Allows splitting of light from intensifier using semitransparent mirror (**beam splitter**)
- 90% of light goes to the film and 10% of light to the TV camera
- aperture used to control amount of light falling on camera faceplate (and to control the degree of vignetting)



# Fiber Optic Coupling

- Advantages:
  - More effective in light collection
  - Less geometrical distortion
- Disadvantage:
  - cannot record image directly from image tube (All recording done from TV)



# Types of TV cameras used

## ➡ **VIDICON TV** camera

- ➡ improvement of contrast
- ➡ improvement of signal to noise ratio
- ➡ high image lag


## ➡ **PLUMBICON TV** camera

- ➡ lower image lag (**follow up of organ motions**)
- ➡ higher quantum noise level
- ➡ **suitable for cardiology**




# How does TV camera operate

- The camera tube is a vacuum tube , that has a photoconductor (anode) in which the light from intensifier output screen is incident.
- In a PLUMBICON tube, this layer is made of lead oxide, whereas in a VIDICON, antimony trisulphide is used
- Electron beam is produced by heating a photocathode , and is focused to the photoconductor by focusing coils so that the electron beam moves across the surface of the TV camera tube in a series of lines
- The surface of the photoconductor scanned with the electron beam results in a flowing current inside the signal plate
- The amount of current is related to the amount of light incident on the conductor.
- The current (video signal) is then transferred to the CRT monitor for display




# CRT Monitor

- The scanning electron beams of the TV camera and the CRT monitor move in synchronism.
- The current of the scanning beam in the TV monitor, is related to that in the TV camera.
- Consequently, the brightness of the image on the TV monitor is proportional to the amount of light falling on the corresponding position on the TV camera.




## *Some points regarding the CRT monitor*


- 1) TV Interlacing: **frame** scanned in two passes, first the odd lines then the even
  - Advantage: avoids flicker
  - each pass takes 1/60th second → One frame takes 1/30th second



- 2) progressive scanning
  - used on newer systems, lines scanned in order
  - no interlacing



- CRT Vertical Resolution: proportional to # of scan lines
- CRT Horizontal Resolution: Depends on modulated signal frequency = cycles/scan line



# Digital fluoroscopy

- Advantages


- 1- computer image storage
- 2-image processing (e.g edge enhancement , geometrical inversion , mapping to look-up tables to ↑contrast)

This processing can be done in real time

# Digital fluoroscopy



## 1) with use of ordinary video camera:

- Signal from video camera → ADC → digital format → displayed on the computer



## 2) With use of CCD camera

- Linked to the intensifier output phosphor by fiber-optic coupling
- CCD camera is formed of Input screen ( amorphous silicon) divided into pixels (1024 x 1024)
- Each pixel acts as a small capacitor
- Charge collected  $\propto$  intensity of light falling on the pixel
- Read-out of pixel charges is so rapid (allow image sequence of 30 frames/s or faster)
- Advantages :
  - produce signal in digital format
  - High dynamic range (12-bit image depth)






# Automatic brightness control (ABC)

- In fluoroscopy manual Kv and mA control is not practical (region of the patient imaged changes rapidly)
- Aim of ABC: automatic adjustment of Kv & mA to keep the output screen light intensity to the required level


Input signal of ABC may be:

- 1-light intensity of the output screen of the image intensifier with optic coupling using a photodiode (older systems)

Area of output screen used is the central area (area of clinical interest), some systems allow operators to choose it manually




2- signal from camera (more in modern systems):  
Signal change feedback to the generator  
→ automatic adjustment of Kv & mA



## Anti-isowatt curves:

- ABC is controlled by different types of programmed curves

**Curve A:** increase in radiological thickness  $\rightarrow \uparrow K_v$  &  $\uparrow mA$   
(maximum power limit = 400 W = 400 J/s)



**Curve B** = image quality weighted curve:

Tube potential is kept between 60-70 Kv , which is optimum for Iodine imaging (k-edge = 33 keV)


So that when the radiological thickness increase  $\rightarrow \uparrow$  mA

At power of 400 W (power rating) tube can't sustain  $\rightarrow$  if more thickness  $\rightarrow \uparrow$  Kv &  $\downarrow$  mA (so that power rating is not exceeded)

Advantages : good image contrast

Disadvantages :  $\downarrow$  kv  $\rightarrow \uparrow$  dose

Uses: contrast studies , angiography



**Curve C:** dose weighted curve:


Increase of radiological thickness  $\rightarrow$   $\uparrow$  Kv rapidly up to maximum Kv that can be sustained by the tube (kv rating)

If more thickness  $\rightarrow$  maximum Kv is used with increasing mA (= high kv curve)

Advantages:  $\uparrow$  Kv  $\rightarrow$  low patient dose

Disadvantage: less image quality

Uses: pediatrics



# Automatic gain control

- Alternative to ABC , but adjust brightness by adjusting the gain of TV system
- Disadvantage : either will  $\uparrow$  image noise or will cause unnecessary dose
- Example:
  - fluoroscopy is used at certain Kv and mA
  - When thickness decrease  $\rightarrow$  kv and mA are stable and AGC is in charge  $\rightarrow$  image brightness is adjusted but the kv and mA are high for the thickness imaged  $\rightarrow$ 
    - 1-unnecessary dose***
    - 2-less noisy image***

*Think what will happen if the thickness increase?*

# Fluoroscopy dose rates

Fluoroscopy dose rate depends on:

- 1- DQE of image intensifier (detection efficiency is high for CsI = 60%)
- 2- level of noise that can be accepted:
  - Input dose rate of image intensifier is usually adjusted to the range of 0.1-0.5  $\mu\text{Gy/s}$  (about 0.25)
  - The lower the input dose , the greater is the noise
  - ABC will give lowest brightness possible at which compensation by TV gain will cause accepted noise



### 3- Magnification (field size) :


magnification  $\rightarrow$   $\downarrow$  of the minification gain  $\rightarrow$  dose is  $\uparrow$  to keep the same brightness  $\rightarrow$   $\uparrow$  ESD

But the area scanned is also decreased  $\rightarrow$  dose area product is unchanged (theoretically)

Yet in practice:


compensation is not only done by  $\uparrow$  dose , but also by  $\uparrow$  TV gain

So that dose area product is reduced to some extent during magnification





#### 4- pulsed fluoroscopy: two types:

- A) primary control pulsed fluoroscopy:
  - X-ray generator switching on and off on regular intervals
- B) Grid controlled tube (direct control):
  - Circular tube electrode between filament and target with –ve voltage (2kv) is switched on and off → electrons pass and then repel and then pass ... etc.
  - Advantages: more precise switching with good control of X-ray



- At rate of 25-30 pulses/s  $\rightarrow$  the eye is not able to detect pulses
- Most of continuous fluoroscopy are actually pulsed
- Each image is retained on the monitor until the following is displayed
- at low pulse rate:
  - 1- no flickering
  - 2- image lag
  - 3- dose rate falls (but not exactly in proportion to the decrease in pulse rate – equipment design dependant)



# Fluoroscopy and ESD

- Entrance surface dose (ESD): can never exceed 100 mGy/min (or 50 mGy/min) for any field of view in fluoroscopy
- In practice with modern equipment : it is in the range of 10-30 mGy/min


# Fluoroscopy recorded images

- Fluoroscopy is usually used for positioning (e.g. biopsy) , it is rarely diagnostic (real-time , and poor image quality)
- For diagnosis, images must be recorded and stored by some way in high quality during fluoroscopy operation

## 1- non digital fluoroscopy

### A) spot films with film screen:

- Cassette holder (with film screen inside) is parked outside the radiation beam
- On pressing the switch to record the image → cassette is driven between the grid and the image intensifier and exposure is made (with delay of about 1 s.)
- Multiple successive exposures can be made in single film using collimator (e.g. barium swallow)





## B) spot films with fluorography

- Definition: recording images produced by image intensifier output screen
- Mirror is introduced between output screen of the image intensifier and the TV camera to take spot films using static camera
- Film size = 105 mm [30% of films produced in method (A)]
- X-ray tube operated in radiographic mode (about 300 mA)
  - i.e. High mA (with short exposure time) → decrease noise

N.B: Cine camera may be used instead of static camera (e.g. cardiac studies)

Using pulsed X-ray exposure (pulse rate = camera frame rate)



## 2- spot images in digital fluoroscopy

- Always done using fluorography

### A- Last image hold = fluorograb:

- Computer stores the last frame of fluoroscopy (continue to appear in the monitor after switching off the X-ray)
- Advantages: operator can inspect the image without additional radiation
- Disadvantages: low quality image

### B- digital spot images:



- Single shot image taken with  $\uparrow$  mA to decrease the noise
- Required image intensifier dose to provide adequate image quality = 0.5 – 5  $\mu$ Gy (exact value depends on magnification degree and clinical requirement)
- If 0.5  $\mu$ Gy is required to produce the spot image  $\rightarrow$  equal to dose produced by 2 seconds of fluoroscopy screening (remember Image intensifier dose rate)
- N.B: Sequential images can be recorded (total number is dependant on computer memory e.g. angiography)

# Spatial resolution of fluoroscopy:

## 1- in image intensifier itself:

- Resolution is limited by spread of light in the output phosphor
- At this step : resolution = 4-5 lp/mm (somewhat better for magnified views)
- Spatial resolution  $\propto$  image size at the input face of image intensifier

i.e. : Resolution is defined in terms of the image size at the input face of image intensifier





## 2- resolution seen in the monitor display:

More limited than step 1

### a) Using CCD :

CCD consists of 1024 x 1024 pixels

So that : for 350 mm FOV  $\rightarrow$  pixel size =  
290  $\mu\text{m}$   $\rightarrow$  spatial resolution = 1.7  
lp/mm


Yet , actually , it is not more than 1.2  
lp/mm)

Improved to 3 lp/m for magnified view

### b) Using Older vacuum tube camera:



Similar resolution as CCD

*What are the factors affecting resolution here?*





# Noise in fluoroscopy

- Significant feature of fluoroscopy (why?)
- Quantum sink corresponds to the photons absorbed in the input screen i.e. nothing will  $\uparrow$  S/N ratio after this step ( e.g. camera gain)
- To decrease the noise:
  - 1) Increase input dose rate by  $\uparrow$  mA
  - 2) frame averaging: adding signals from successive frames
    - Advantages:
      - a-  $\downarrow$  noise (equivalent to increasing exposure time) with higher ability to see low contrast details
      - B- smoothing the image
    - Disadvantages: blurring if high level of motion between frames
    - $\uparrow$  number of frames added  $\rightarrow \downarrow$  image noise  $\rightarrow \uparrow$  blurring



Www.ScoreTraining.net



Frame 1

Frame 1+2

Frame 1+2+3

Frame 1+2+3+4


# Veiling glare



- Causes: scattering effect in the image intensifier:
  - 1- light scattering in the output window (mainly)
  - 2- X-ray and light scatter in the input phosphor
  - 3- electron scattering in the tube itself
- The larger is the image intensifier , the greater is the Veiling glare
- Effect:
  - 1- decrease the contrast because dark regions of the image appear lighter
  - 2- makes the central area of the image appear brighter than the periphery = vignetting

N.B: other causes of vignetting :

- A) periphery of the image is displayed over larger area of the input screen
- B) poor periphery focusing



# Geometrical distortion

- Types:

- 1- pin cushion distortion:

Magnification towards the edges of the image



Due to curvature of the input screen

- 2- S- type distortion:

Impose curvature on the straight line features


Due to ambient magnetic fields influencing path of electrons from input to output screens

*Both are insignificant when imaging complex shapes within the body*





# Testing the spatial resolution

- Using grid test
- Test is carried out without any attenuating materials in the beam
- ↓ kv is used to minimize the scatter
- Detect any deterioration in the image intensifier focusing





# Testing the contrast resolution

- under standard fluoroscopy conditions (70 kv , 1-2 mm copper filter)
- Using low contrast test object (e.g. Leeds test object) which is a flat disc , 6mm in thickness , 200-300 mm in diameter
- Contain groups of circular inserts of  $\uparrow$  Z,
- Example: 12 groups of circular details
- Each group has diameter A , and Each group is composed of 9 disks
- Thickness of the disks of each group increase progressively to produce varying levels of contrast
- Observer is required to count the number of details that can be seen in each group
- Usually it is not possible to see more than 6 details in any group , and more than 1 detail in the smallest diameter group



- $H_t(A) = [C_t(A) \times A^{\frac{1}{2}}]^{-1}$
- $H_t(A)$  : threshold detection index of group A
- $C_t(A)$  : threshold contrast of group A (minimum visible contrast)
- A: diameter of the group
- $H_t(A)$  is plotted against A to detect imaging performance (by comparing it to the previous tests)



# 3) Flat Panel Detectors

- Flat panel detectors are becoming more and more used in fluoroscopy.
- Flat panel detectors will replace the image intensifier-TV based fluoroscopy systems
- Indirect conversion (with an intermediate step) most common





# The use of flat panel detectors instead of image intensifier in fluoroscopy

	Flat panel detectors	Image intensifier with CCD camera
DQE	60-65 %	60-65% (both use CsI)
Dynamic range and contrast	Very high (14 bit depth) → better contrast	CCD has good contrast resolution (12 bit depth) , yet , contrast is limited due to veiling glare → contrast ratio is not better than 30:1
Pixel size	150 $\mu\text{m}$ (resolution = 3lp/mm)	400 $\mu\text{m}$ (resolution = 1-1.2 lp/mm)
magnification	Does not affect spatial resolution (pixel size is the same)	Resolution improve with magnification (for the smallest field of view resolution of both systems are comparable
Geometrical distortion	-ve	+ve
CT scanner	Can be used in multislice CT scanners	Can be used but limited due to circular FOV and geometrical distortion
Dr.Yossef Gamal	Www.ScoreTraining.net	

# Types of fluoroscopy tables

1) TUBE ABOVE THE TABLE = *Over couch Tubes*



## 2) TUBE UNDER THE TABLE = under *couch* X-ray tube



### 3) C-ARM FLUOROSCOPY:

Tube and intensifier move together around the patient



# Digital subtraction angiography

- Objectives:

- 1) produce image of contrast filled vessels in isolation from other tissues

This will increase vessels clarity and decrease dose of contrast

- 2) rapid image succession to avoid patient movement





**actual mask image**



**actual angiographic image**




**↓ B - A (Computer process)**



**angiographic sub-traced image**


- Process:
  - A) two non contrast images are taken (first : to stabilize x-ray factors , second: used as mask image (image A))
  - B) post-contrast image is taken (image B)
  - C) A is subtracted from B on a pixel to pixel basis → result in subtraction image C or D (shows the filled vessels only)
  - D) sequence of subtracted images can be recoded , based on the initial mask (and may be viewed in real time)
  - E) signals in the contrast and mask images are first converted to their logarithms before subtraction , and then are converted back into intensity values for display





## ***Problem(1): motion misregistration:***

- *Cause:* movement between frames
- Particularly at boundaries between high contrast objects (e.g. at bone edges)
- *Solution:* pixel shifting:
  - retrospective adjustment of the mask image
  - This can be done only over the full area of the image (no differential movement within the field of view)
  - Can be done manually or using automated techniques




## ***Problem (2): increased noise***


- *Cause* : subtraction cause that the number of pixels contributing to each pixel is decreased
- *Solution*: DSA requires increased mA



**2  $\mu$ Gy per frame**



**15  $\mu$ Gy per frame**



**24  $\mu$ Gy per frame**

***Problem (3): length of anatomy imaged is usually greater than image intensifier field of view (e.g. lower limb angio)***

- Solution (1):
  - several mask images are acquired
  - separate contrast administration for each region
- Solution (2):
  - Several mask images are acquired
  - Single contrast image is acquired along the full leg length with longitudinal movement of the table to track the progress of contrast
  - Appropriate mask image for each table position is used in the subtraction
  - Advantage: ↓amount of contrast medium used
  - Disadvantages: increase likelihood of movement between mask and contrast images






# Rotational angiography

- 90° rotation x-ray tube and intensifier about the patient → mask images are taken at several angles
- Post-contrast images are acquired at the same angles
- DSA are produced at each angle
- Advantages:

1-identification of optimum angle to view the vessel

2- three dimensional reconstruction



# Dual energy subtraction

# Dual energy subtraction

1) image A is taken at low Kv  $\rightarrow$  high contrast between bone and soft tissue

2) image B is taken at high Kv  $\rightarrow$  contrast is reduced

3) Subtraction process:

- $B - A \rightarrow$  minimize visualization of bones and improves soft tissue contrast (e.g. chest x-ray to remove rib densities)
- $A - B \rightarrow$  displays bony structures in greater details




Standard image



Soft tissue image




Bone image



# Imaging Screen MTF

- improves with magnification
- Why?
  - magnification of an object of given frequency reduces frequency seen by screen




# Computed tomography





# Radiography vs. CT


- Radiography
  - Structures are superimposed on film (Viewed through underlying / overlying structures)
  - Thin highly-attenuating objects appear to be same density as thicker low-attenuating object.



- CT
  - Cross-sectional image  
i.e. Image computed from beam intensity measurements through only slice of interest
  - View anatomy without looking through underlying / overlying structures → improves the contrast
  - Uses tightly collimated beam → minimizes scattered radiation, improves contrast



# Pixels and Voxels




# Pixels & Voxels

- Pixel is 2D component of an image



- Voxel is 3D cube of anatomy




- Voxel depth

  - same as slice thickness

- CT reconstruction calculates attenuation coefficients of Voxels

- CT displays these attenuation coefficients as gray shades of Pixels



# CT Number = Hounsfield Units


Determines the grade of gray which is displayed in each pixel  
Depends on the attenuation coefficient of the voxel

$$\text{CT \#} = 1000 \times \frac{(\mu_t - \mu_w)}{\mu_w}$$

Where:

$\mu_t$  = linear attenuation coefficient for tissue in voxel

$\mu_w$  = linear attenuation coefficient for water



**Example 1: voxel contains water ( $\mu_t = \mu_w$ ):**

**Example:  $\text{CT\#} = 1000 \times (\mu_t - \mu_w) / \mu_w$  voxel contains air ( $\mu_t \approx 0$ ):**

$$\text{CT\#} = 1000 \times (0 - \mu_w) / \mu_w = 1000 \times (-1) = -1000$$

- Notes:

1- air and water are used for calibration of CT number scale of the scanner

2- CT number depends on  $\mu$  of the tissue

i.e. depends on : Z, Density , kv used & filtration


3- range of CT number values can be displayed is :  
-1024 to + 3071

Explanation:

Pixel depth of the CT scanners usually = 12


i.e. levels of gray displayed =  $2^{12} = 4096$

Tissue	CT Number (HU)
<b>Bone</b>	<b>+1000</b>
Liver	40-60
White mater	-20 to -30
Grey mater	-37 to -45
Blood	40
Muscle	10-40
Kidney	30
CSF	15
<b>Water</b>	<b>0</b>
Fat	-50 to -100
Air	-1000






# Volume averaging effect

- Cause: The CT number stored in each voxel represents average attenuation coefficient in the voxel
- Results:
  - High contrast object can be seen even if it is smaller than the pixel (small low contrast details will not be visualized)
  - thin **high** contrast structure will be visualized larger than normal e.g. oblique enhanced vessel





Depends on: slice thickness : the thinner the slice, the less is the volume averaging




# Display of CT (Window and level)


- Gray shade assigned to each pixel value (CT #)
- If the image displayed as such: image with very bad contrast (very high latitude)

47

93




- Windowing:
  - Assignment of densities to pixel values according to clinical interest
  - does not disturb original pixel values in memory
  - Automatically or manually controllable
    - window
    - level



Window: 400  
Level: 500

# CT # Window Change

Window: 800  
Level: 500





Dr.Yossef Gamal

Www.ScoreTraining.net

All pixels with HU <301 will appear black

All pixels with HU >900 will appear

# CT Window Change




- *The narrower is window width to the greater is the image contrast and the greater is the ability to reveal low contrast structures*
- *The narrower is window width to the less is the displayable latitude*

# CT equipment

# CT components


- gantry :
  - Provides structural framework for all other components
- Couch: Moving longitudinally
- x-ray tube : Produces x-rays
- Filters and Collimators
- Detectors



## CT planes:

X and y = trans-axial plane

Z = longitudinal plane = plane of table movement








# **CT scanners generations**




# • First Generation Scanners

- Formed of X-ray source and a single detector
- Tube produce a pencil beam
- Translation/Rotation type:
  - Translation: both tube and detector move across the scanning plane to acquire series of measurements
  - Rotation: then tube and detector is rotated by  $1^\circ$  and process is repeated
- 5 minutes to gather enough information for one slice
- Tube was only able to rotate 180 degrees





# Second Generation

- Fan-shaped x-ray beam
- Bank of detectors (30 or more), still insufficient to cover full cross section of the patient
- Rotation translation movement is still needed
- 20 seconds per slice
- 180 degree rotation





# Third Generation

- Fan-shaped x-ray beam
- Large number of small detectors arranged in an arc to cover the full cross section of the patient → linear translation is not required
- Complete 360 degree rotation of tube and detectors
- Rotate/Rotate movement
- One rotation = one slice
- Time reduced to 1 sec per slice
- This type will be discussed in more details (Helical and multislice scanners are based on its idea)



# Fourth Generation

- Rotate – stationary scanner
- Rotating x-ray tube
- Detectors are arranged stationary ring completely surrounding the patient
- Advantages:
  - More detectors stability
  - Simpler reconstruction
  - Outer part of fan beam always passes outside the patient → during each rotation every detector is able to measure the unattenuated beam → this measurement is used to adjust calibration throughout the scan
- Disadvantages:
  - ↑ number of detectors (6 folds)
  - ↑ dose due to ↑ distance between patient and detectors





# Fifth Generation = electron beam scanners

- Electron beam focused on target ring which cover  $210^\circ$  arc around the patient
- $\rightarrow$  X-ray is produced & collimated and detected by detectors ring above the patient forming  $216^\circ$  ring
- Advantage:
  - No rotation of mechanical parts  $\rightarrow$  electron beam is swept across the full arc in less than 50 ms. (slice time)
- Used in cardiology

